
**Operating Instructions | Betriebsanleitung | Mode d'emploi | Manuale d'uso |
Manuale d'uso | Instrucciones de manejo | Инструкция по эксплуатации |
操作说明 | 取扱説明書 |**

YDK03MS | YDK04MS

**Density Determination Kit | Dichtebestimmungsset | Dispositif de détermination de la
masse volumique | Dispositivo per la determinazione della densità | Kit densitométrico |
Комплект для определения плотности | 密度測定工具 | 密度測定キット**

English – page 3

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Contents

1 About this Document	4	5 Operation	13
1.1 Validity	4	5.1 Method used to Determine Density	13
1.2 Other Applicable Documents	4	5.2 Sources of Error and Options for Correction	14
1.3 Symbols Used	4	5.3 Density Determination	17
1.3.1 Warning Information in Operation Descriptions	4	5.3.1 Determining Density of Solid Objects	17
1.3.2 Other Symbols	4	5.3.2 Determining the Density of Solid Objects with a Density of Less Than 1 g/cm ³	18
2 Safety	5	5.3.3 Determining the Density of Liquids	21
2.1 Intended Use	5	5.4 Density Values	22
2.2 Glass Components	5		
3 Device Description	6	6 Formulas	23
3.1 YDK03MS	6		
3.2 YDK04MS	7		
4 Installation	8		
4.1 Scope of Delivery	8		
4.2 Unpacking the Device	8		
4.3 Installing the Product	8		
4.3.1 Attaching the YDK03MS to the Balance	8		
4.3.2 Attaching the YDK04MS to the Balance	10		
4.3.3 Selecting the Beaker and Immersion Device	11		
4.3.4 Installing the Glass Plummet	12		
4.3.5 Installing the Thermometer	12		

1 About this Document

1.1 Validity

These instructions are part of the product. These instructions apply to the product in the following versions:

Product	Type
Density determination kit for Cubis® laboratory balances with a readability of 0.1 mg and 0.01 mg	YDK03MS
Density determination kit for Cubis® laboratory balances with a readability of 1 mg	YDK04MS

1.2 Other Applicable Documents

- ▶ In addition to these instructions, please observe the following documents:
 - Instructions for the laboratory balance used

1.3 Symbols Used

1.3.1 Warning Information in Operation Descriptions

NOTICE

Denotes a danger with risk that property damage may result if the risk is **not** avoided.

1.3.2 Other Symbols

- ▶ Required action: Describes actions which must be carried out.
- ▷ Result: Describes the result of the actions carried out.

2 Safety

2.1 Intended Use

This product is used to determine the density of solid objects or liquids. For the determination of the density of solid objects, a test liquid must be used.

This product must be operated with Sartorius laboratory balances with a suitable readability.

If the laboratory balance has a density determination program: Density values can be calculated automatically using the density determination kit and density determination program.

If the density determination program is used to determine density: The density must be determined as described in the laboratory balance operating instructions.

The product is intended exclusively for use in accordance with these instructions. Any further use beyond this is considered **improper**.

2.2 Glass Components

The glass plummet, beakers, and thermometer are made of glass. Broken glass can cause injuries.

- ▶ If the glass plummet, beakers, or thermometer become damaged: Replace the glass plummet, beakers, or thermometer.

3 Device Description

3.1 YDK03MS


Fig. 1: YDK03MS components

No.	Name	No.	Name
1	Beakers (Ø 90 mm and Ø 55 mm)	5	Glass plummet
2	Frame	6	Sieve
3	Thermometer retainer clip, installed upon delivery	7	Sample holder
4	Thermometer	8	Bridge

3.2 YDK04MS


Fig. 2: YDK04MS components

No.	Name	No.	Name
1	Beakers (Ø 90 mm and Ø 55 mm)	5	Glass plummet
2	Frame	6	Sieve
3	Thermometer retainer clip, installed upon delivery	7	Sample holder
4	Thermometer	8	Bridge

4 Installation

4.1 Scope of Delivery

Item	Quantity
Beaker (Ø 90 mm and Ø 55 mm)	2
Frame	1
Thermometer retainer clip	1
Thermometer	1
Glass plummet	1
Sieve	1
Sample holder	1
Bridge	1

4.2 Unpacking the Device

Procedure

NOTICE

Unpacking can damage components!

► This product consists of delicate components. Take care when unpacking all parts.


4.3 Installing the Product

4.3.1 Attaching the YDK03MS to the Balance


Prerequisites

A Cubis® laboratory balance with a readability of 0.1 mg or 0.01 mg.

Procedure

- ▶ Remove the following parts from the balance:
 - Weighing pan (1)
 - Pan support (2)

- ▶ Insert the frame into the weighing chamber. The wedge-shaped opening at the top of the frame must face in the direction from which the sample holder (sieve/glass plummet) will be placed into the frame.


- ▶ Place the bridge which will hold the beaker through the frame on the balance.

4.3.2 Attaching the YDK04MS to the Balance


Prerequisites

A Cubis® laboratory balance with a readability of 1 mg.


Procedure

- ▶ Remove the weighing pan (1) from the balance.

- ▶ Insert the frame into the weighing chamber. The wedge-shaped opening at the top of the frame must face in the direction from which the sample holder (sieve/glass plummet) will be placed into the frame.


- ▶ Place the bridge which will hold the beaker through the frame on the balance.

4.3.3 Selecting the Beaker and Immersion Device


Procedure

- ▶ Select the beaker or immersion device.
The selection of the beaker or immersion device depends on the sample which is to be determined.


Overview of the beaker and immersion device:

- For determining the density of solid objects whose density is higher than the immersion liquid: A beaker with \varnothing 90 mm and a sample holder is required

- For determining the density of solid objects whose density is lower than the immersion liquid: A beaker with \varnothing 90 mm and a sample pan is required


- For determining the density of liquids: A beaker with \varnothing 55 mm and a glass plummet is required

4.3.4 Installing the Glass Plummet

A wire is attached to the glass plummet by a wire eyelet. The wire eyelet must be suspended from the hook on the glass plummet.

Procedure


- ▶ **NOTICE** The wire can break if it is bent! Pull the glass plummet out of the packaging by the glass eyelet. Do **not** bend the wire on the glass plummet.
- ▶ Loop the wire eyelet of the glass plummet over the hook of the retainer.
- ▶ Attach the hook of the glass plummet to the frame.

4.3.5 Installing the Thermometer

Procedure

- ▶ Place the thermometer in the beaker.
- ▶ If the thermometer needs to be attached: Attach the thermometer to the rim of the beaker using the retainer clip.

5 Operation

5.1 Method used to Determine Density

"Archimedes' Principle" is used to determine the density of a solid object with this measuring device:

- An object immersed in liquid will be subject to an upward buoyant force. This force is equal to the weight of the liquid displaced by the object.
- Using a hydrostatic balance, which enables you to weigh the solid object in air and in water, it is possible to determine the density of a solid object if the density of the buoyancy medium is known:

$$\rho = \frac{W(a) \cdot \rho(fl)}{W(a) - W(fl)}$$

or

Determine the density of a liquid if the volume of the immersed object is known.

$$\rho(fl) = \frac{G}{V}$$

Where:

ρ	density of the solid object
$\rho(fl)$	density of the liquid
$W(a)$	weight of the solid object in air
$W(fl)$	weight of the solid object in liquid
G	buoyancy of the immersed object
V	volume of the immersed object

5.2 Sources of Error and Options for Correction

NOTICE

Incorrect procedure when adjusting the balance!

- ▶ Remove the sample holder of the density determination kit prior to making adjustments.
-

The formula used to determine the density of solid objects is sufficient to obtain an accuracy of one to two decimal places. Depending on the accuracy required, the following error and correction factors will need to be considered:

- effect of temperature on the density of the buoyancy liquid
- air buoyancy when weighing in air
- changes in the immersion depth of the iron pan when immersing the sample
- adhesion of the liquid to the span wire of the iron pan
- air bubbles sticking to the sample

Some of the errors can be corrected mathematically. It is therefore necessary to measure the temperature of the liquid and correct the liquid density accordingly and to define the inner diameter of the vessel holding the liquid.

Effect of Temperature on Liquid Density

The density of the buoyancy liquid is temperature dependent. The change in density per °C temperature change is as follows:

- 0.02% for distilled water
- 0.1% for alcohols and hydrocarbons

The change can show up in the third decimal place during the density determination.

To correct the liquid density in terms of temperature, proceed as follows:

- measure the temperature of the liquid using the thermometer supplied.
- The density of water can be found in the table in Chapter "5.4 Density Values," page 22 and can be used for ρ (fl).

Air Buoyancy

Depending on the temperature, humidity and air pressure, a 1 cm³ volume of air will have a weight of around 1.2 mg. When weighing in air, the object experiences a corresponding buoyancy per cm³ of its volume. The error that results if the air buoyancy is not allowed for shows up in the third decimal place and should therefore be corrected.

The buoyancy force is taken into account in the following formula:

$$\rho = \frac{W(a) \cdot [\rho(fl) - \rho(a)]}{W(a) - W(fl)} + \rho(a)$$

Where $\rho(a) = 0.0012 \text{ g/cm}^3$ = density of air under normal conditions (temperature 20°C, pressure 101.325 kPa).

Immersion Depth

The pan for holding and/or immersing the sample during weighing in liquid is attached rigidly to two wires and is immersed approximately 30 mm below the surface of the liquid. Since the balance is tared prior to each measurement, the additional buoyancy from the submerged part of the measurement setup is not factored into the determination of the density.

When weighing in liquid, a volume of liquid corresponding to the volume of the sample body gets displaced. This causes the fastening wires on the pan to be immersed deeper and to generate additional buoyancy, creating an error in the density determination.

The following formula will correct the error:

$$\rho = \frac{W(a) \cdot [\rho(fl) - \rho(a)]}{0.99983 [W(a) - W(fl)]} + \rho(a)$$

As the correction factor is determined exclusively by the geometry of the device, it must be ensured that only the supplied vessel with the larger diameter (90 mm) is used to determine the density of a solid object. Chapter "6 Formulas," page 23 shows how this correction factor is derived.

Adhesion of the Liquid to the Wire

When immersing the sample holder (or sieve) in the buoyancy liquid, liquid travels up the wire because of adhesion forces and generates a few additional milligrams of weight.

Since the sample holder (or sieve) is in the buoyancy medium both when weighing in air and when weighing in liquid, and the balance is tared at the beginning of each measurement, the influence of the liquid meniscus can be ignored.

To reduce the surface tension and friction of the liquid on the wire, add three drops of a surfactant (Mirasol Antistatic or a conventional detergent) to the distilled water in the vessel.

Because of the buoyancy liquid travelling up the wire, the weight value may continue to slowly change after the "g" has appeared. For this reason, the weight value should be read as soon as the "g" appears.

Air Bubbles

The measurement errors which occur as a result of air bubbles sticking to the sample can be evaluated as follows. If the air bubble has a diameter of 0.5 mm, this will produce an additional buoyancy of less than 0.1 mg when weighing in water. If the air bubble has a diameter of 1 mm, the additional buoyancy will be around 0.5 mg, and if the diameter is around 2 mm, roughly 4.2 mg. It is therefore imperative that larger air bubbles are taken off with a fine brush or similar.

Moisture can also be added in advance in a separate container.

5.3 Density Determination

5.3.1 Determining Density of Solid Objects

Preparation

In the following description, distilled water is used.

Procedure

- ▶ Place the beaker with the larger diameter (\varnothing 90 mm) on the center of the bridge
- ▶ Fill the beaker until the distilled water is approximately 5 mm below the rim.
- ▶ Add three drops of surfactant to the distilled water.
- ▶ Attach the thermometer to the rim of the beaker using the retainer clip.
- ▶ Clean the sample holder with a solvent (paying particular attention to the immersed wires) and hook on to the frame.

Measurement Procedure

The sample weight must be determined in the air.

Procedure

- ▶ Tare the balance.
- ▶ Place the sample on the frame weighing pan and weigh it.
- ▶ Make a note of the weight value $W(a)$.

Determining the Buoyancy

$$G = W(a) - W(fl)$$

Procedure

- ▶ Tare the balance with the sample on the frame weighing pan.
- ▶ Place the sample in the sample holder. Remove the iron pan from the measuring device and ensure that **no** additional air bubbles become attached upon re-immersion in the liquid; If possible, add the sample directly using tweezers or similar.
- ▶ Make a note of the absolute value of buoyancy G , with a minus sign in front.

Calculating the Density

Procedure

- ▶ Read off the temperature.
- ▶ Locate the density value ρ (fl) in the table using the temperature you have read off.
- ▶ Calculate the density using the following formula:

$$\rho = \frac{W(a) \cdot [\rho(\text{fl}) - 0.0012 \text{ g/cm}^3]}{0.99983 G} + 0.0012 \text{ g/cm}^3$$

$W(a)$ and G in g; $\rho(\text{fl})$ in g/cm^3 $G = W(a) - W(\text{fl})$

5.3.2 Determining the Density of Solid Objects with a Density of Less Than 1 g/cm^3

There are two different methods for determining the density of solid objects with a density of less than 1 g/cm^3 .

Method 1

Distilled water is still used as the buoyancy liquid. However, the iron pan is replaced by the inverted sieve pan (sieve).

To determine the buoyancy of the sample, first place it on the surface of the water and then immerse it using the sieve which was previously removed.

It is also possible to use tweezers or similar to place the sample directly under the sieve pan (without removing the sieve from the frame).

If the buoyancy of the substance to be measured is greater than the weight of the sieve, the weight of the sieve must be increased by adding additional weight to the frame weighing pan.

Method 2

A liquid with a density which is lower than that of the solid object which is to be measured must be used as a buoyancy medium. Good results have been achieved with ethanol (up to a density of approx. 0.8 g/cm³).

The negative effect of the liquid's surface tension on the results is less noticeable when using ethanol than when using distilled water. It is therefore not necessary to add surfactants.

Valid safety precautions must be taken when working with ethanol.

The second method should be used if the density of the solid object only differs slightly from that of distilled water. As the sample is suspended in water, measuring errors may occur if the first method is used. It also makes sense to use the second method to determine the density of a granulated substance. In this case, it would be difficult to place the entire sample under the sieve when using the first method.

Do not use ethanol if the sample could be attacked or dissolved by it.

Preparation

In the following description, distilled water is used.

Procedure

- ▶ Place the beaker with the larger diameter (\varnothing 90 mm) on the center of the bridge
- ▶ Fill the beaker until the distilled water is approximately 5 mm below the rim.
- ▶ Add three drops of surfactant to the distilled water.
- ▶ Attach the thermometer to the rim of the beaker using the retainer clip.
- ▶ Clean the sieve with a solvent (paying particular attention to the immersed wires) and hook on to the frame.

Measurement Procedure

The sample weight must be determined in the air.

Procedure

- ▶ Tare the balance.
- ▶ Place the sample on the frame weighing pan and weigh it.
- ▶ Make a note of the weight value $W(a)$.

Determining the Buoyancy

$$G = W(a) - W(fl)$$

Procedure

- ▶ Tare the balance with the sample on the frame weighing pan.
- ▶ Place the sample in the sample holder. Remove the iron pan from the measuring device and ensure that **no** additional air bubbles become attached upon re-immersion in the liquid; If possible, add the sample directly using tweezers or similar.
- ▶ Make a note of the absolute value of buoyancy G , with a minus sign in front.

Calculating the Density

Procedure

- ▶ Read off the temperature.
- ▶ Locate the density value $\rho(fl)$ in the table in Chapter "5.4 Density Values," page 22 using the temperature you have read off.
- ▶ Calculate the density using the following formula:

$$\rho = \frac{W(a) \cdot \rho(fl)}{0.99983 G} + 0.0012 \text{ g/cm}^3$$

$W(a)$ and G in g; $\rho(fl)$ in g/cm^3 $G = W(a) - W(fl)$

5.3.3 Determining the Density of Liquids

Preparation

Procedure

- ▶ Place the beaker with the smaller diameter (\varnothing 55 mm) on the center of the bridge.
- ▶ Attach the thermometer to the rim of the beaker using the retainer clip.

Measurement Procedure

Procedure

- ▶ Suspend the disk with the glass plummet (hanging on one wire) from the frame.
- ▶ Tare the balance.
- ▶ Fill the beaker with the liquid to be tested until it is 10 mm above the glass plummet.

Determining the Buoyancy $G = W (a) - W (fl)$

The negative weight displayed by the balance corresponds to the buoyancy acting on the glass plummet in the liquid.

Procedure

- ▶ Make a note of the buoyancy G , with a minus sign in front.
- ▶ Read off and make a note of the temperature.

Calculating the Density

Procedure

► Calculate the density using the following formula:

$$\rho \text{ (fl)} = \frac{G}{V}$$

G in g; V in cm³

The glass plummet included in the density determination kit has a volume of 10 cm³.

It is easy to obtain the actual density of the liquid (in g/cm³); Shift the decimal point in the balance display one place the left.

5.4 Density Values

Density of H₂O (Gas-Free Water) at Temperature T (in °C) in g/cm³

T/°C	0.0	0.1	0.2	0.3	0.4	0.5	0.6	0.7	0.8	0.9
10	0.99970	0.99969	0.99968	0.99967	0.99966	0.99965	0.99964	0.99963	0.99962	0.99961
11	0.99960	0.99959	0.99958	0.99957	0.99956	0.99955	0.99954	0.99953	0.99952	0.99951
12	0.99950	0.99949	0.99947	0.99946	0.99945	0.99944	0.99943	0.99941	0.99940	0.99939
13	0.99938	0.99936	0.99935	0.99934	0.99933	0.99931	0.99930	0.99929	0.99927	0.99926
14	0.99924	0.99923	0.99922	0.99920	0.99919	0.99917	0.99916	0.99914	0.99913	0.99911
15	0.99910	0.99908	0.99907	0.99905	0.99904	0.99902	0.99901	0.99899	0.99897	0.99896
16	0.99894	0.99893	0.99891	0.99889	0.99888	0.99886	0.99884	0.99883	0.99881	0.99879
17	0.99877	0.99876	0.99874	0.99872	0.99870	0.99869	0.99867	0.99865	0.99863	0.99861
18	0.99860	0.99858	0.99856	0.99854	0.99852	0.99850	0.99848	0.99846	0.99844	0.99842
19	0.99840	0.99838	0.99837	0.99835	0.99833	0.99831	0.99829	0.99826	0.99824	0.99822
20	0.99820	0.99818	0.99816	0.99814	0.99812	0.99810	0.99808	0.99806	0.99803	0.99801
21	0.99799	0.99797	0.99795	0.99793	0.99790	0.99788	0.99786	0.99784	0.99781	0.99779
22	0.99777	0.99775	0.99772	0.99770	0.99768	0.99765	0.99763	0.99761	0.99758	0.99756
23	0.99754	0.99751	0.99749	0.99747	0.99744	0.99742	0.99739	0.99737	0.99734	0.99732
24	0.99729	0.99727	0.99725	0.99722	0.99720	0.99717	0.99714	0.99712	0.99709	0.99707
25	0.99704	0.99702	0.99699	0.99697	0.99694	0.99691	0.99689	0.99686	0.99683	0.99681
26	0.99678	0.99676	0.99673	0.99670	0.99667	0.99665	0.99662	0.99659	0.99657	0.99654
27	0.99651	0.99648	0.99646	0.99643	0.99640	0.99637	0.99634	0.99632	0.99629	0.99626
28	0.99623	0.99620	0.99617	0.99615	0.99612	0.99609	0.99606	0.99603	0.99600	0.99597
29	0.99594	0.99591	0.99588	0.99585	0.99582	0.99579	0.99576	0.99574	0.99571	0.99568
30	0.99564	0.99561	0.99558	0.99555	0.99552	0.99549	0.99546	0.99543	0.99540	0.99537

References:

- NIST ITS-90 Density of Water Formulation for Volumetric Standards Calibration, Volume 97, Number 3, Chapter 3.3
Density of Air-Free Water ρ (kg m⁻³) = 999.85308 + 6.32693 x 10⁻² t - 8.523829 x 10⁻³ t² + 6.943248 x 10⁻⁵ t³ - 3.821216 x 10⁻⁷ t⁴
- DAKS-DKD-LF-8-2 calibration liquid water – DaksS, annex 1, table A

6 Formulas

To gain a better understanding of the process, the theory behind the formulas and the correction factor is explained here.

Basic Principles

$$\text{Density} = \frac{\text{Mass (g)}}{\text{Volume (cm}^3\text{)}}$$

Archimedes' principle states:

An object immersed in liquid will be subject to a buoyant force (G). This force is equal to the weight of the liquid displaced by the object. The volume of the immersed object V (s) is equal to the volume of the displaced liquid V (fl).

The following are determined:

1. Weight in air W (a)
2. Buoyancy of the object in the liquid (G)

The density of an object is:

$$\rho = \frac{\text{mass of object}}{\text{volume of object}} = \frac{W (a)}{V (s)} = \frac{W (a)}{V (fl)}$$

If the density ρ (fl) of the displaced liquid is known:

$$V (fl) = \frac{\text{Mass (fl)}}{\rho} = \frac{G}{\rho (fl)}$$

Thus:

$$\rho = \frac{W (a) \cdot \rho (fl)}{G}$$

Calculation

The density of a solid object is calculated based on the ratio:

$$\rho : W (a) = \rho (fl) : [W (a) - W (fl)]$$

Thus:

$$\rho = \frac{W (a) \cdot \rho (fl)}{W (a) - W (fl)}$$

$$W (a) - W (fl) = G = \text{buoyancy of the sample}$$

The density of a liquid is determined from the buoyancy of the immersed object which has a defined volume.

$$\rho (fl) = \frac{G}{V}$$

Where:

ρ	density of the solid object
$\rho (fl)$	density of the liquid
$W (a)$	weight of the solid object in air
$W (fl)$	weight of the solid object in liquid
G	buoyancy of the immersed object
V	volume of the immersed object

Corrections

The following must be allowed for when correcting the density determination:

- The air buoyancy experienced by the sample when weighed in air.
Where $\rho (a) = 0.0012 \text{ g/cm}^3$ = density of air under normal conditions (temperature 20°C, pressure 101.325 kPa); which results in:

$$\rho = \frac{W (a) \cdot [\rho (fl) - \rho (a)]}{W (a) - W (fl)} + \rho (a)$$

- The immersion of the wires of the sample holder or sieve

When using this specific density determination kit, the buoyancy $G = [W(a) - W(fl)]$ must be multiplied by the factor 0.99983 (corr). Expanded formula:

$$\rho = \frac{W(a) \cdot [\rho(fl) - \rho(a)]}{[W(a) - W(fl)] \cdot \text{corr}} + \rho(a)$$

This factor allows for the buoyancy of the wires which are submerged deeper when the sample is introduced.

How this Correction Factor Is Derived:

The buoyancy caused by the submerged wires depends on the height "h" by which the liquid rises when the sample is immersed.

Here, the sample volume $V(pr)$ corresponds to the liquid volume $V(fl)$.

The sample volume is determined by measuring the buoyancy. Hence:

$$V(pr) = V(fl)$$

or

$$\frac{W(a) - W(fl)}{\rho(fl)} = \frac{\pi \cdot h \cdot D^2}{4}$$

$$\text{Then } h = \frac{4 \cdot [W(a) - W(fl)]}{\rho(fl) \cdot \pi \cdot D^2}$$

The buoyancy "A" caused by the immersed wires is:

$$A = 2 \cdot \frac{\pi \cdot d^2}{4} \cdot h \cdot \rho(fl)$$

When "h" is used:

$$\rho = \frac{2 \cdot \pi \cdot d^2 \cdot 4 \cdot [W(a) - W(fl)] \cdot \rho(fl)}{4 \cdot \rho(fl) \cdot \pi \cdot D^2}$$

$$A = 2 \cdot \frac{d^2}{D^2} \cdot [W(a) - W(fl)]$$

To allow for the buoyancy of the wires, subtract the buoyancy "A" caused by the immersed wires from the buoyancy determined for the sample:

$G = W(a) - W(fl)$. The buoyancy value "A (corr)" which is to be used in the calculation is then: $G - "A"$.

$$A(\text{corr}) = [W(a) - W(fl)] - 2 \cdot \frac{d^2}{D^2} \cdot [W(a) - W(fl)]$$

$$A(\text{corr}) = \left[1 - 2 \cdot \frac{d^2}{D^2} \right] \cdot [W(a) - W(fl)]$$

In this density determination kit, the beaker with the larger diameter (\varnothing 90 mm) and an immersion device with two wires with a diameter of 0.7 mm are used to determine the density of solid objects.

When the values $d = 0.7$ mm and $D = 90$ mm are used, the correction factor is:

$$1 - 2 \cdot \frac{0.7^2}{90^2} = \mathbf{0.99983}$$

When using devices with other dimensions, the correction factor must be recalculated accordingly.

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Phone: +49.551.308.0
www.sartorius.com

The information and figures contained in these instructions correspond to the version date specified below.

Sartorius reserves the right to make changes to the technology, features, specifications and design of the equipment without notice. Masculine or feminine forms are used to facilitate legibility in these instructions and always simultaneously denote the other gender as well.

Copyright notice:

This instruction manual, including all of its components, is protected by copyright. Any use beyond the limits of the copyright law is not permitted without our approval. This applies in particular to reprinting, translation and editing irrespective of the type of media used.

© Sartorius Germany

Last updated:
10 | 2017

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Sartorius Lab Instruments GmbH & Co. KG