
**Installation Instructions | Installationsanleitung | Mode d'emploi |
Istruzioni per l'uso | Instrucciones de funcionamiento**

Sartorius YRD03Z

**Remote Display | Zweitanzeige | Afficheur supplémentaire |
Display supplementare | Display adicional**

English page 3

In cases involving questions of interpretation, the German-language version shall prevail.

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Contents

	Page
Intended Use	3
Warning and Safety Information	3
Equipment Supplied	3
Installation	4

Intended Use

YRD03Z is a remote display which can be directly connected to Sartorius scale (see next page). The remote display is equipped with a back-lit 7-digit LCD. A number of symbols can be displayed in addition to the “g” (gram) weight unit symbol, including: kg, ct, lb, oz, ozt, tl, GN, dwt, mg, m, k, t, b, m, pcs and %. Unit symbols that cannot be displayed are indicated by an “o”.

The remote display is operated through an RS-232 interface. The cable supplied has standard 25-pin connectors and is either approx. 23 cm or 93 cm long (see package contents).

Note:

- △ The remote display may not be used for legal metrology unless certain measures are carried out. Please contact your Sartorius representative for details.

Symbols


The following symbols are used in these instructions:

- indicates operating instructions
- indicates for operating instructions that should only be carried out under specific circumstances
- > describes what occurs after an operation is carried out
- indicates items in a list
- △ indicates a hazard

Warning and Safety Information

- △ The remote display must not be connected to scales which are operated with batteries.
- △ When cleaning, do not permit liquids to enter the remote display.
- △ Do not use the mains connector socket on the remote display. The power is supplied via the data interface of the scale.

Equipment Supplied

- 1× Remote Display YRD03Z
- 1× Connection cable, length: 93 cm.

- Connection cable for installation, length: 23 cm:
- 1× Connection cable, length: 23 cm.
- 1× Fixing plate
- 4× Countersunk screws
- 4× Threaded bushes

Installation


Installation of the Remote Display (Cable: 93cm) on Sartorius Models of the Following Series:

AC, BJ*, BL*, BP, CP, CPA, ED, GC, GK, GM*, GP, LA,
LC, LE, LP, MC, ME, RC, SC, SE, TE,


Note:

The scale must be equipped with a 25-pin data interface (RS-232).


- * Possibility of retrofitting a data interface as an option (see also the operating instructions for the scale). For connection to other series, please contact your Sartorius representative.

- Disconnect the scale from the power supply:
Unplug power cord.

- Unplug the mains unit jack plug.


Installing the Data Cable (93cm)

- Remove the cable (93cm) from the package supplied.
- Plug the 25-pin plug into the interface of the scale and secure it.


- Plug the 9-pin plug into the interface of the scale and secure it.

Note:

Select the appropriate menu setting
(see respective scale operating instructions).

- Reconnect the scale to the power supply.


Installation of the Remote Display (Cable: 23cm) on Sartorius Models of the Following Series:

ED2201/-CW, ED2202S/-CW, ED3202S/-CW,
ED4202S/-CW, ED5201/-CW, ED6202S/-CW,
ED8201/-CW, GW 3202, GW 6202, GW 7201.

Note:


The scale must be equipped with a 25-pin data interface (RS-232).

ED and GW models with serial number 22505142 or higher. For connection to other series, please contact your Sartorius representative.


- Disconnect the scale from the power supply:
Unplug power cord.

- Unplug the mains unit jack plug.

- Remove the cable (23 cm) from the package supplied.


- Remove the weighing pan.

- Turn over the scale and the remote display and place on a soft surface (e.g.: foam rubber mat).

- Carefully press the threaded bushings into the base plate of the scale and into the remote display with your thumbs.


- Plug the data cable (23 cm) into the data output of the scale and into the remote display (9-pin) and secure by tightening the screws.


- Reconnect the mains unit jack plug.

Note:


- ⚠ Do not use the mains connector socket on the remote display.
The power is supplied to the remote display via the data cable.
- ⚠ For technical reasons, the scale must remain voltage-free in this position. Do not plug in the mains unit. Do not apply voltage to the scale!


- Place the fixing plate (the countersunk holes point upwards) so that the center foot of the scale is inserted through the hole.
- Push the remote display onto the scale.
- Insert the countersunk screws into the threaded bushings and tighten them in order to fasten the fixing plate.

- Turn over the scale and the remote display.
- > The remote display is now installed.

- Bring the scale with the remote display to the place where it is to be set up.
- Reconnect the power supply. Plug in the mains unit.
- Switch on the scale (ON/OFF switch).
- Calibrate/adjust the scale before use (see operating instructions).

- Alternatively, the remote display may be mounted on a normal photographic tripod (not supplied).

Adapting Data Output for Operation of the Remote Display

The data output of the scale must be set to “automatically without stability”. Check the setting of the scale and change it if necessary.

Settings:

See the appropriate operating instructions, Section: Configuration, scale operating menu/scale operating program or the interface description for the scale.

Menu change and settings according to model, starting from the default settings:

Series:	Code
AC, BP, CP, CPA, GC, GP, LC, LE, LP, MC, SC, RC,	6 1 4
BJ, BL, GM, TE	6 1 3

For the following series, select the menu path:

Series:

ED, GK and GW

Path:

```

Setup  Report
      (Print)
          Print
          (manual/automatic)
          Automatic without stability
  
```

Series:

LA

Note:

The particular menu setting (automatic weight display) must be taken into account for each individual application.

Path:

Setup
(print
output)

Data Output

to interfaces

Serial Communication
(PERIPHERALS)

Automatic
Display Value Output

Note:

The particular menu settings (automatic display value output) must be taken into account in each individual application.

Series:

ME, SE

Path:

Setup

Print Output

Output to Interfaces

Serial Communication
(PERIPHERALS)

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Copyright by Sartorius,
Goettingen, Germany.
All rights reserved. No part
of this publication may
be reprinted or translated in
any form or by any means
without the prior written
permission of Sartorius.
The status of the information,
specifications and illustrations
in this manual is indicated
by the date given below.
Sartorius reserves the
right to make changes to
the technology, features,
specifications and design of the
equipment without notice.

Status:
October 2011,
Sartorius Weighing Technology GmbH,
Goettingen, Germany